

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
CHT100110	Treatment Chair - Greiner Multiline Next IT V2 - Navy Blue Vinyl - 220kg	\$ 8,249.00	\$ 2,474.70	1
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	2
CHP198970	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Large - 160kg ^	\$ 3,689.00	\$ 1,844.50	2
SCT670540	Scooter - Pride - Celebrity Deluxe XL - 4 Wheel Red (with Rear bag + Cane Holder)	\$ 6,379.00	\$ 3,189.50	1
MWS442600	Wheelchair Folding Push/Pull - 440mm Leverdrive - Left Hand	\$ 4,629.00	\$ 2,314.50	1
CHP197706	Chair Powerlift Recline - Aspire Signature - Space Saver - Small - Heritage - 145kg	\$ 2,099.00	\$ 1,049.50	2
CHP199500	Chair Powerlift Recline - Air Comfort - Airwing - Lge 600x510x710mm - 1 Motor	\$ 3,829.00	\$ 1,148.70	1
CHP198975	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Large - 160kg	\$ 3,689.00	\$ 1,844.50	1
LSS390592	Standing & Transport Aid - Shifty Transfer Assist	\$ 3,099.00	\$ 929.70	1
LIC360760	FreeSpan Rail 5m - with Likorall Trolley - 200kg	\$ 3,049.00	\$ 1,524.50	1
LSS393360	RollOn - Max 160 kg (350 lbs).	\$ 2,800.00	\$ 1,400.00	1
PTB535402	SuperShape Back Support - 15in Wide x 20in High with adjustable shape and contour system	\$ 2,739.00	\$ 1,369.50	1
WAT907201	Avalo AX Articulative Arm with Pole Mount (needs 1x A&E Mount 12149K)	\$ 447.00	\$ 134.10	6
BEB045850	Bed - Aidacare AC3 V2 - Single - 3 Function - Teak Bed Ends	\$ 2,450.00	\$ 1,225.00	1
CHP197741	Chair Powerlift Recline - Aspire Signature - Dual Action - Small - Heritage Vinyl - 155kg	\$ 2,269.00	\$ 1,134.50	1

Item Number	Item Description	WAS	NOW	QTY
BEA007480	Alarm - Cura - Bed Monitor Plug Pack 12V DC (for 3500 Active Care Consol) (replacement)	\$ 29.10	\$ 8.73	1
BEA007516	Nurse Call Test Box	\$ 120.00	\$ 36.00	1
BEA007520	Alarm - Monitor - Mono - To Suit Nurse Call Systems	\$ 184.00	\$ 55.20	1
BEA014605	Bed Cradle - Aspire - Fixed - Treated Steel ^	\$ 69.00	\$ 20.70	1
BEA029705	Self Help Pole - Aidacare AC/Ward Beds	\$ 184.00	\$ 92.00	1
BEB045850	Bed - Aidacare AC3 V2 - Single - 3 Function - Teak Bed Ends	\$ 2,450.00	\$ 1,225.00	1
BEP054300	Pillow Cover - Pyramid Ruffled Slip only	\$ 54.00	\$ 27.00	1
BTC066300	Bedpan and Bowl Carrier - 4HC (& KCare) - for Shower commodes	\$ 210.00	\$ 63.00	2
BTC067108S	Shower Commode Elevating Legrest - Aspire - LEFT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC067109S	Shower Commode Elevating Legrest - Aspire - RIGHT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC067110S	Shower Commode Elevating Legrest - Aspire - LEFT - Tilt in Space Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC067115S	Shower Commode Elevating Legrest - Aspire - RIGHT - Tilt in Space Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC069530	Shower Commode Seat - 4HC - Closed Front w/ Pan Carrier Runners - 550mm Wide	\$ 469.00	\$ 140.70	1
BTC078300	Anti Tip Bars - for K Care Shower Commodes - OLD - Compatible w Candy Cane Leg	\$ 165.60	\$ 49.68	1
BTC080360	Foot Board - K Care - Shower Commode - Grey Plastic	\$ 124.00	\$ 37.20	2
BTC087636	Anti Tippers - Invacare - for Aquatec Ocean & VIP Shower Commodes	\$ 123.00	\$ 36.90	2
BTC089180	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 500mm - Closed Front	\$ 461.00	\$ 138.30	1
BTR101751	Rail 1500mm Tube 32mm diameter BGE	\$ 73.80	\$ 22.14	4
BTR101771	One T-piece, one wall flange with washer & screw BGE	\$ 19.10	\$ 5.73	4
BTR101780	One 30 deg bend, one wall flange with fixings	\$ 31.40	\$ 9.42	7
BTR101801	One 90 deg bend, one wall flange with fixings BGE	\$ 31.40	\$ 9.42	1
BTR101830	4 way fitting, wall flange includes fixings	\$ 31.10	\$ 9.33	11
BTR101851	Rail 900mm Tube 32mm diameter BGE	\$ 43.70	\$ 13.11	1
BTR101865	Beige End Set for 1 Corner: Two elbows, one std wall flange, 1corner flange	\$ 24.00	\$ 7.20	2
BTR101881	Toilet Rail Holder with fixing for grabrail or wall	\$ 23.00	\$ 6.90	8
BTR101885	Internal Corner fitting, one corner wall flange BGE	\$ 35.50	\$ 10.65	3
BTR101890	Gravity toggles Large or small with centreing sleeve	\$ 8.50	\$ 2.55	16
BTR101962	Large End Ring	\$ 14.10	\$ 4.23	10
BTR101963	Connector Ring	\$ 17.60	\$ 5.28	25
BTR103000	Rail 1500mm - Ripple Finish - Beige	\$ 36.10	\$ 10.83	5
BTR103150	Ramp Grab Rail - Modular - Tubing 1000mm - Brushed S/Steel	\$ 48.60	\$ 14.58	1
BTR103320	Toilet roll holder - Stainless Steel for 4H11L and 4H11R	\$ 18.70	\$ 5.61	1
BTR103410	Grab Rail - Modular - Elbow 90 degree horizontal LH - Brushed S/Steel	\$ 65.00	\$ 19.50	7
BTS105400	Bath Board Alton Bariatric 190kg	\$ 285.00	\$ 85.50	1
BTS106000	Shower Stool w/ Arms & Padded Seat - K Care - 550mm Seat - Zinc	\$ 778.40	\$ 233.52	2
BTS150100	Shower Chair - Freedom - HBA407 (Oval Tube) - Aluminium HD Shower Chair - 540mm - 250kg	\$ 275.00	\$ 82.50	2
BTS151110	Shower Chair - Saljol - Spa - Medium - Dark Grey - 150kg	\$ 220.00	\$ 66.00	1
BTS151120	Shower Chair - Saljol - Spa - Small - Dark Grey - 150kg	\$ 220.00	\$ 66.00	1
BTT149700	Toilet Support Rail Toilet Roll Holder - for Devon Rail	\$ 74.00	\$ 22.20	4
BTT156500	Over Toilet Aid Toilet Seat Lid - for U Frame OTA	\$ 23.70	\$ 7.11	2
BTT156501	Over Toilet Aid Toilet Seat - Replacement for OTA U Frame	\$ 23.70	\$ 7.11	1
BTT163110	Over Toilet Aid - Freedom - HBA370-HD - Aluminium HD Dignity Split Seat Chair - 560mm - 200kg	\$ 325.00	\$ 162.50	3
CHP197055	Chair Powerlift Recline - Pride - C101 Riser Recliner Pecan - 510x510x700mm - 1 Motor	\$ 1,391.50	\$ 695.75	1
CHP197706	Chair Powerlift Recline - Aspire Signature - Space Saver - Small - Heritage - 145kg	\$ 2,099.00	\$ 1,049.50	2
CHP197741	Chair Powerlift Recline - Aspire Signature - Dual Action - Small - Heritage Vinyl - 155kg	\$ 2,269.00	\$ 1,134.50	1
CHP197771	Head/Armrest Cover Set - Aspire Signature Recliners - Lagoon	\$ 46.30	\$ 13.89	1
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	2
CHP198970	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Large - 160kg ^	\$ 3,689.00	\$ 1,844.50	2
CHP198975	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Large - 160kg	\$ 3,689.00	\$ 1,844.50	1
CHP198991	Advanced Immersion Backrest Cushion - Medium - Navy - Aspire COVE Chair	\$ 427.00	\$ 128.10	1
CHP199500	Chair Powerlift Recline - Air Comfort - Airwing - Lge 600x510x710mm - 1 Motor	\$ 3,829.00	\$ 1,148.70	1
CHP213300	Hilite Chair Rear Adj Legs w/ Wheels - Vanilla/Champ	\$ 201.00	\$ 60.30	1
CHP216300	Chair Footrest - Princess Chair/Bed & Deluxe - Air Comfort	\$ 477.00	\$ 143.10	1
CHP216353	Chair Tray - Air Comfort - Deluxe V2 - Large	\$ 385.00	\$ 115.50	1
CHP217300	Chair Tray Large - Tilt / Deluxe Bed/Daily Chair - Air Comfort	\$ 262.00	\$ 78.60	3
CHP227100	Chair Powerlift Recline - Aspire Verona - Mushroom - Single Motor - 130kg ^	\$ 1,659.00	\$ 829.50	1

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
CHT100110	Treatment Chair - Greiner Multiline Next IT V2 - Navy Blue Vinyl - 220kg	\$ 8,249.00	\$ 2,474.70	1
CON240650	Cello Form Pads AMD Maxi Small (Pkt 20)	\$ 57.00	\$ 17.10	9
CON241000	Cello AMD Form Anatomical Shaped Pads PE (Maxi Abs) Liliac (Pkt 20) (2450ml)	\$ 40.00	\$ 12.00	7
CON243300	TENA Pads Comfort Maxi (Pkt 28)	\$ 53.00	\$ 15.90	5
CON249960	TENA Pads Flex Maxi Extra Large (Pkt 21)	\$ 107.00	\$ 32.10	5
DCS100180	Dementia Care - Orientation Signs - Kitchen - Foam Board, Red	\$ 71.30	\$ 21.39	1
DCS100240	Dementia Care - Orientation Signs - Toilet - Foam Board, Yellow	\$ 71.30	\$ 21.39	1
DIN265830	Tumbler Autoplas 230ml Natural	\$ 5.49	\$ 1.65	37
DLG272215	Lighting - Outdoor Sensor Light	\$ 71.00	\$ 21.30	1
DLK276110	Cutlery - Kings Double lightweight foam handle w strap	\$ 41.20	\$ 12.36	1
FBS280560	Elastic Rib Belt - Body Assist - Mens - Large	\$ 68.00	\$ 20.40	1
FBS283980	Wrist/Hand - Support - BodyAssist - Slip-On Elastic - Small Beige	\$ 34.90	\$ 10.47	1
FBS284650	Elbow Sleeve - BodyAssist - Slip-On Elastic - Bge - Med	\$ 24.90	\$ 7.47	1
FBS301159	Hand - Glove - Isotoner - Therapeutic pressure gloves - Open Finger - Sml	\$ 63.00	\$ 18.90	1
FBS301531	Knee Hinged - Oppo - ROM - Elastic Wrap ROM - O/Patella - Sml	\$ 112.00	\$ 33.60	1
FP308454	Exercise - Theraband - 152cm Black Special Heavy	\$ 14.10	\$ 4.23	13
GMD300320	Sphygmomanometer - Digital Semi Auto Hand Held	\$ 198.00	\$ 59.40	1
GMD300375	Sphygmomanometer - Aneroid Wall Model 6in Aneroid Dial	\$ 317.00	\$ 95.10	1
GSS340251	Cushion Cover - Waterproof cover for Latex Multi-use Roll	\$ 58.00	\$ 17.40	1
LAC353030	Linen Bag Waterproof Black	\$ 70.00	\$ 21.00	1
LIC353360	Ceiling Hoist - Prism - Bracket End Stop - Long w/Charger	\$ 416.00	\$ 124.80	1
LIC353370	Ceiling Hoist - Prism - Charging Dock for Transactive Xtra	\$ 374.00	\$ 112.20	1
LIC356010	CP200 (Prism) Bottom Cover	\$ 43.70	\$ 13.11	3
LIC359150	Carriage Hook for Multirall or Likorall	\$ 187.00	\$ 56.10	1
LIC359160	Adapter for Multirall or Likorall	\$ 43.70	\$ 13.11	9
LIC359580	Hang-Up Hand Control Hanger (10 pcs) (standard on Likorall 242/243/250 and Multirall 200)	\$ 46.20	\$ 13.86	1
LIC360670	HandControl HB80 - Liko	\$ 374.00	\$ 112.20	1
LIC360760	FreeSpan Rail 5m - with Likorall Trolley - 200kg	\$ 3,049.00	\$ 1,524.50	1
LSS371635	Lifter Sling - Oxford Quick Fit Deluxe Padded Head Support Sling Small	\$ 404.00	\$ 202.00	1
LSS372700	Lifter Sling - Access Toileting - Kerry - Large	\$ 288.00	\$ 144.00	1
LSS373310	Lifter Sling - Invacare Hammock - Large	\$ 224.00	\$ 112.00	1
LSS373400	Sling General Purpose Fab Small	\$ 327.00	\$ 163.50	1
LSS374100	Sling General Purp Fab XXL (S.M.S.)	\$ 285.00	\$ 142.50	1
LSS374912	Sling General PurposeXXLarge & Head Support	\$ 311.00	\$ 155.50	1
LSS374916	Sling - Hygiene with Head Support Polyester - Small	\$ 327.00	\$ 163.50	2
LSS376100	Sling - Supreme Hygiene - Fabric Small	\$ 327.00	\$ 163.50	2
LSS376422	Sling - Hygenic - Fabric - Large	\$ 327.00	\$ 163.50	1
LSS377754	Lifter Sling - Allegro - Care Head Support Mesh - Small	\$ 402.00	\$ 201.00	1
LSS378800	Lifter Sling - Standing Lifter - Kerry - Large	\$ 297.00	\$ 148.50	1
LSS379000	Lifter Sling - Seating Sling - Kerry - Large	\$ 342.00	\$ 171.00	1
LSS379900	Lifter Sling - Pivot Frame Polyester w/ Keyhole - Kerry - Small	\$ 364.00	\$ 182.00	1
LSS380700	Lifter Sling - Gen Purp Mesh - Kerry - Large LR406NE3	\$ 297.00	\$ 148.50	1
LSS380900	Lifter Sling - Gen Purp Polyester - Kerry - Small	\$ 689.00	\$ 344.50	1
LSS382000	Lifter Sling - Gen Purp Mesh & Head Support - Kerry - Large LR407NE3	\$ 349.00	\$ 174.50	1
LSS382100	Lifter Sling - Gen Purp Mesh & Head Support - Kerry - Extra Large LR407NE4	\$ 669.00	\$ 334.50	1
LSS382800	Lifter Sling - Full Body Mesh w/ Div Leg - Kerry - Large LR408NE3	\$ 376.00	\$ 188.00	1
LSS390592	Standing & Transport Aid - Shifty Transfer Assist	\$ 3,099.00	\$ 929.70	1
LSS393360	RollOn - Max 160 kg (350 lbs).	\$ 2,800.00	\$ 1,400.00	1
LSS394500	Sling - Liko HygieneSling w/ Belt - Poly - Reinforced leg supports - Lge - 200kg	\$ 519.00	\$ 259.50	1
LSS394750	Sling - Solo Highback Sling - Liko - Med - 2x10 pcs - No Wash - 200kg	\$ 1,629.00	\$ 814.50	1
MWS442600	Wheelchair Folding Push/Pull - 440mm Leverdrive - Left Hand	\$ 4,629.00	\$ 2,314.50	1
MWS449205	Wheelchair Folding Light - Action 3NG - Self Propelled - 16in (405mm) - 24in solid tyres	\$ 1,329.00	\$ 664.50	1
PTA493750	HipSaver Slim Fit Veterans Kit Female X-Large	\$ 529.00	\$ 370.30	1
PTA494005	HipSaver Slim Fit Veterans Kit Male Small	\$ 529.00	\$ 370.30	1
PTA504300	Hip Protector Pants Female Size 10 (S.M.S.)	\$ 34.00	\$ 23.80	2
PTA504700	Hip Protector Pants Female Size 18 (S.M.S.)	\$ 34.00	\$ 23.80	3
PTA504900	Hip Protector Pants Female Size 22 (S.M.S.)	\$ 34.00	\$ 23.80	3
PTA505400	Hip Protector Pants Male Small (S.M.S.)	\$ 37.00	\$ 25.90	3
PTA508480	WChair Armpad Trough Liner - Pride - Large - Gel	\$ 203.00	\$ 142.10	2
PTA511900	HipSaver NH Veterans Kit Male Small	\$ 529.00	\$ 370.30	1
PTA512438	HipSaver - Quick Change - Veterans Kit 2 (6-Pants, 2-Pads) - Female - Med	\$ 529.00	\$ 370.30	1
PTA512800	HipSaver - QuickChange - Starter Kit (3-Pants, 1-Pads) - Female - Med	\$ 269.00	\$ 188.30	2
PTA535040	Tortoise Turning and Positioning System - Molnycke - Standard	\$ 999.00	\$ 699.30	1
PTB533882	WChair Back - JAY J3 - Mounting Bracket - Short Left	\$ 34.40	\$ 17.20	3
PTB533886	WChair Back - JAY J3 - Extra Mounting Hardware	\$ 358.00	\$ 179.00	2
PTB534152	Spex - Vigour Mid backrest upgrade to contouring cell-cube system	\$ 789.00	\$ 394.50	1
PTB535402	SuperShape Back Support - 15in Wide x 20in High with adjustable shape and contour system	\$ 2,739.00	\$ 1,369.50	1
PTB535457	Spex Back Support Quad Mount Hardware set for second wheelchair	\$ 519.00	\$ 259.50	1
PTB538850	WChair Back - Dreamline Contour Backrest - Fixed - 18in x 22in	\$ 1,449.00	\$ 724.50	1
PTB539520	WChair Back - Dreamline Support Backrest - Quick Release - 16in x 18in	\$ 1,529.00	\$ 764.50	1
PTB539870	WChair Back - Dreamline Support S3 Backrest - Quick Release - 18in x 22in	\$ 1,529.00	\$ 764.50	1
PTB540580	WChair Back - Dreamline BAC System - 18in x 22in	\$ 659.00	\$ 329.50	1
PTC551572	Cushion Ride Forward 16in x 16in Complete	\$ 769.00	\$ 384.50	1
PTC561274	Cushion - Jay X2 - 17 x17in	\$ 759.00	\$ 379.50	1
PTC566701	Cushion - J2 Deep Contour Cushion 24x18in	\$ 1,549.00	\$ 774.50	1

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
PTC571902	Cushion - JAY - Lite 14x14in	\$ 509.00	\$ 254.50	1
PTC571970	Cushion - JAY - Lite 15x16in	\$ 509.00	\$ 254.50	1
PTC576500	Cushion Cover - JAY J2 - 15.5x16in - Recline Air Exch	\$ 126.00	\$ 37.80	2
PTC594920	Cushion - ROHO Hybrid Elite 1R 16x18in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594935	Cushion - ROHO Hybrid Elite 1R 18x16in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC596500	Cushion Cover - ROHO Standard HP 18.5x22in	\$ 153.00	\$ 45.90	1
PTC600600	Cushion Synergy Structure 18 x 18"	\$ 772.00	\$ 386.00	1
PTC608015	Cushion - Promedicare Domino - 10 x 12in w/ 3D cover	\$ 1,619.00	\$ 809.50	1
PTC620900	Cushion - Dreamline Contour Cushion - 16in x 18in	\$ 979.00	\$ 489.50	1
PTC622035	Cushion - Axiom - Standard - 20 x 18	\$ 157.00	\$ 78.50	2
PTM636345	ROHO Mattress Enclosure Cover	\$ 434.00	\$ 303.80	2
PTM636550	Mattress Overlay - Suzrik Bubble Mat	\$ 257.00	\$ 128.50	1
PTM636625	AirLo - Thermal Hyperstretch Microfibre Coverlet	\$ 494.00	\$ 148.20	1
PTM638020	Mattress - Aspire HYBRID ACUTE - Single - 9 Cell - 1980x895x170mm - 230kg (Pump Sep: PTM638030)	\$ 1,939.00	\$ 969.50	1
PTS640880	Stockings - JOBST UltraSheer - Knee - 30-40mmHg - C/Toe - Med - Blk	\$ 102.00	\$ 30.60	2
PTS641310	Stockings - JOBST Relief - Thigh w/ Grip Top - 20-30mmHg - Unisex - O/Toe - Med	\$ 130.00	\$ 39.00	1
PTS641974	Stockings - Venosan 6001 (Nyl/Lyc) - Knee - 18-22mmHg - Unisex - O/Toe - Short - Sml - Blk	\$ 107.00	\$ 32.10	1
SCA668311	Scooter Anti Rattle Clamp for "L" Piece	\$ 17.10	\$ 5.13	1
SCA668415	Scooter Acc - Invacare - Colour Kit Green	\$ 168.00	\$ 50.40	1
SCA668420	Scooter Acc - Invacare - Colour Kit Blue	\$ 168.00	\$ 50.40	1
SCT670540	Scooter - Pride - Celebrity DeluxeXL - 4 Wheel Red (with Rear bag + Cane Holder)	\$ 6,379.00	\$ 3,189.50	1
TSE672527	Transfer Tube	\$ 36.10	\$ 10.83	1
TSE674240	Defibrillator Table - Aspire Emergency Stretchers	\$ 323.00	\$ 161.50	2
WAA686250	Crutch Acc - Aspire - Grey Padded Hand Grip	\$ 20.70	\$ 6.21	2
WAA695615	Crutch & Cane Holder - Aspire VOGUE INDOOR Walkers	\$ 12.50	\$ 3.75	2
WAF750150	Seat Walker - Freedom - Stride - Red - 135kg	\$ 369.00	\$ 184.50	5
WAF751100	Seat Walker - Saljol - Page Indoor Walker - Small - Charcoal - 150kg	\$ 380.00	\$ 114.00	1
WAS748110	Walking Stick - Freedom - Cobra Cane - 710-940mm - 110kg	\$ 56.00	\$ 16.80	1
WAS748130	Walking Stick - Freedom - Crook Handle - Aluminium w/Black Handle - 770-970mm - 110kg	\$ 26.40	\$ 7.92	1
WAS748170	Walking Stick - Freedom - Crook Handle - Bronze w/Black Handle - 770-970mm - 110kg	\$ 23.70	\$ 7.11	1
WAS748180	Walking Stick - Freedom - Crook Handle - Floral - 770-970mm - 110kg	\$ 28.50	\$ 8.55	1
WAS748250	Walking Stick - Freedom - Offset Handle - Black - 770-970mm - 110kg	\$ 24.00	\$ 7.20	2
WAS748260	Walking Stick - Freedom - Offset Handle - Bronze - 770-970mm - 110kg	\$ 24.00	\$ 7.20	1
WAS748270	Walking Stick - Freedom - Offset Handle - Chrome - 770-970mm - 110kg	\$ 24.00	\$ 7.20	6
WAT818500	FlexiCart Options - Waste Bin (No Bracket)	\$ 105.00	\$ 31.50	1
WAT820100	IV Pole PH - Kerry	\$ 155.00	\$ 46.50	1
WAT907200	Avalo A&E Mount Bracket Assembly (K) - Universal L&R (single)	\$ 98.00	\$ 29.40	3
WAT907201	Avalo AX Articulating Arm with Pole Mount (needs 1x A&E Mount 12149K)	\$ 447.00	\$ 134.10	6
WCA848170	WChair Elevating Legrest - Advanced - Right - 320-410mm - Aspire Assist2/Transit2/Evoke2/Rehab RS&R	\$ 183.00	\$ 54.90	1
WCA848330	WChair Contoured Padded Overlay - 350mm - Aspire Evoke2	\$ 99.00	\$ 29.70	1
WCA848530	WChair Depth Extension Kit - 50mm - Black - Aspire Rehab RS/Rehab RX	\$ 25.80	\$ 7.74	2
WCA848550	WChair Therapist Kit for Hand Adjustment - Aspire Rehab RS/Rehab RX	\$ 38.30	\$ 11.49	2
WCA858005	WChair Elevating L/rests - alu f/plates - Heavy Duty WChair (20-24in)	\$ 152.00	\$ 45.60	2
WCA858008	WChair Stump Support - Right - Heavy Duty WChair (20-24in)	\$ 340.00	\$ 102.00	1
WCA858350	Lap Belt - Aspire EVOKE Wheelchair	\$ 20.00	\$ 6.00	1
WCA858390	WChair Elevating Legrest - Aspire EVOKE 400/450/500/HD - Black Right	\$ 120.00	\$ 36.00	3
WCA858400	WChair Elevating Legrest - Aspire ASSIST/TRANSIT - Black Left	\$ 110.00	\$ 33.00	1
WCA858462	WChair Anti Tip Bar - Aspire ASSIST/TRANSIT	\$ 25.00	\$ 7.50	1
WCA858610	WChair Belt w/ Plastic Clip Continuous	\$ 93.00	\$ 27.90	3
WCA858823	Headrest Clamp Heavy Duty (New Style) - Black	\$ 268.00	\$ 80.40	1
WCA858920	Side Thoracic Support Clamp 3/4" (19mm) Black	\$ 70.00	\$ 21.00	4
WCA859170	Modular Hand Pad Horn Small Left	\$ 125.00	\$ 37.50	1
WCA859275	WChair Foot Plate Protectors	\$ 62.00	\$ 18.60	1
WCA859416	Start Anti-Tipper Right - Pendulum Style	\$ 167.00	\$ 50.10	2
WCA859423	Stump Support - Left - 20in - Ottobock Start 4 Wheelchairs	\$ 402.00	\$ 120.60	1
WCA859439	Stump Support - Ottobock - For M2S MWC - Left - 16-18in	\$ 323.00	\$ 96.90	1
WCA859441	Stump Support - Ottobock - For M2S MWC - Right - 16-18in	\$ 323.00	\$ 96.90	1
WCA860386	BodyPoint Padded Hip Belt Rear Pull Four Point Compact Push Button Medium	\$ 334.00	\$ 100.20	1